



DEUTSCH

Toilettensitzerhöhung
Gebrauchsanleitung

ENGLISH

Raised toilet seat
Operating instructions

FRANÇAIS

Rehausse WC
Mode d'emploi

ITALIANO

Rialzo per WC
Istruzioni per l'uso

ESPAÑOL

Alzador del asiento del inodoro
Manual de instrucciones

NEDERLANDS

Toiletbrilverhoging
Gebruiksaanwijzing

SVENSKA

Toalettstolsförhöjare
Bruksanvisning

SUOMI

WC-istuimen korotus
Käyttöohje

NORSK

Toalettseteforhøyelse
Bruksanvisning

DANSK

Toiletforhøjelse
Betjeningsvejledning

Ελληνικά

Ανυψωτικό τουαλέτας
Οδηγίες Χρήσεως



Hersteller:


eureha GmbH
Zeppelinstraße 10
D-88353 Kisslegg, Germany
E-Mail info@eureha.de
Internet www.eureha.de




Sicherheitshinweise
Safety instructions
Consignes de sécurité
Indicazioni di sicurezza

Indicaciones de seguridad
Veiligheidsinstructies
Säkerhetsanvisningar

Turvaohjeet
Sikkerhetskhenvisninger
Sikkerhedshenvisninger
Οδηγίες Ασφαλείας



Max. Belastung
Max. Load
Charge max.
Portata max.
Carga máx.
Max. belasting
Max. belastning
Maks. kuormitus
Maks. belastning
Maks. belastning
Μέγιστο Βάρος



Bei Defekt nicht verwenden
Do not use if defective
Ne pas utiliser en cas de défaut
Non impiegarlo se difettoso
No utilizar en males condiciones
Indien defect niet gebruiken
Använd inte produkten om den är trasig
Älä käyttää viallista tuotetta
Må ikke brukes når den er defekt
Må ikke anvendes, hvis den er defekt
Μην το χρησιμοποιείται αν είναι ελαττωματικό





Absturzgefahr
Danger of falling
Risque de chute
Pericolo di caduta
Peligro de caída
Gevaar voor vallen
Risk att ramla
Putoamisvaara
Fare for å falle
Fare for at falde ned
Κίνδυνος Πτώσης

Lieferumfang
Scope of delivery
Volume de livraison
Fornitura

Volumen de entrega
Omvang van de levering
Leveransomfattning




Toimituskokonaisuus
Leveringsomfang
Indeholdt i leveringene
Σκοπός Παράδοσης



Montieren
Mounting
Montage
Montaggio

Montar
Monteren
Montering




Asennus
Montering
Montering
Βάση



Reinigen
Cleaning
Nettoyage
Pulizia

Limpieza
Reinigen
Rengöra

Puhdistaminen
Rengjøring
Rengjøring
Καθαρισμός



Reinigungstemperatur max. 85 °C
Max. cleaning temperature 85 °C
Température de nettoyage 85 °C
Temperatura max. durante la pulizia 85 °C
Temperatura máx. de limpieza 85 °C
Reinigungstemperatur max. 85 °C
Rengjøringstemperatur max. 85 °C
Puhdistuslämpötila maks. 85 °C
Temperatur ved rengjøring maks. 85 °C
Rengjøringstemperatur maks. 85 °C
Μεγ. Θερμοκρασία καθαρισμού 85 °C

Ohne Scheuermittel reinigen
Clean without scouring agent
Ne pas utiliser de produits récurrents pour le nettoyage
Non pulirlo con prodotti abrasivi
No utilizar productos abrasivos para su limpieza
Zonder schuurmiddelen reinigen
Använd inte skurmedel
Puhdista ilman hankausaineita
Bruk ikke skuremiddel ved rengjøring
Rengør uden skuremiddel
Μη χρησιμοποιείτε σύρμα κουζίνας για καθαρισμό

Maschinenwaschbar
Washable in an autoclave
Lavable dans un autoclave
Lavabile in un'autoclave
Lavable en una autoclave
Wasbaar in een autoclaaf
Tvättbar i en autoklav
Pesunkestävä keskelle eräs autoklaavi
Vaskbar i en autoklave
Vaskebjørn i en autoklave
Πλένεται σε συσκευή αποστείρωσης